

## Supply List - Canon PGI-220/225 & CLI-221/226 Cartridges

Website: [www.rjettek.com](http://www.rjettek.com)

Phone #: 800-688-0129 (Ext. 114)

Note: Enter Part# in this search block on website →




Cartridge #	Size	Part No.	Description	Units	Min. Qty Purchase
<b>Ink</b>					
PGI-220/225 Bk	Std	BKP44	OCP Ink Canon BCI-3/23.. PG-240/240XL Pig Blk	oz	4
CLI-221/226 Bk	Std	BK124	OCP Ink Canon CLI-221 Black	oz	4
CLI-221/226 C	Std	C144	OCP Ink Canon CLI-221 Cyan	oz	4
CLI-221/226 M	Std	M144	OCP Ink Canon CLI-221 Magenta	oz	4
CLI-221/226 Y	Std	Y144	OCP Ink Canon CLI-221 Yellow	oz	4
CLI-221/226 Gy	Std	BK123	OCP Ink Canon CLI-221 Gray	oz	4
<b>Chips</b>					
All	Std	1110	Ball Plug Remover	each	1
PGI-220 Bk	Std	0249	Chip, Canon PGI-220 Black, 10pk	10/pk	1 pk
CLI-221 Bk	Std	0244	Chip, Canon CLI-221 Black, 10pk	10/pk	1 pk
CLI-221 C	Std	0245	Chip, Canon CLI-221 Cyan, 10pk	10/pk	1 pk
CLI-221 M	Std	0246	Chip, Canon CLI-221 Magenta, 10pk	10/pk	1 pk
CLI-221 Y	Std	0247	Chip, Canon CLI-221 Yellow, 10pk	10/pk	1 pk
CLI-221 Gy	Std	0248	Chip, Canon CLI-221 Grey, 10pk	10/pk	1 pk
<b>Parts and Tools for Refilling</b>					
All	Std	1904S	Needle, Sharp, 1.5" x 18 gauge	each	1
All	Std	1914	Ball Plug, HP 564, 920, 932, 950, and Many More	each	25
All	Std	1916	Plug, HP 564/920, Many Canon Tanks	each	1
All	Std	1951	Syringe Plastic Luer Lock 30cc Non-Sterile	each	1
All	Std	1953	Syringe Plastic Luer Lock 20cc Non-Sterile	each	1
All	Std	3205-USB	Canon Chip Resetter 5/8/220/221/225/226	each	1
<b>Packaging</b>					
CLI-221/226	Std	1813-OR	Cap, Orange, HP & Canon Small Ink Ports	each	1
PGI-220/225 BK	Std	1814-OR	Cap, Orange, HP & Canon Large Ink Ports	each	1
All	Std	1911	Plug Cover For Various HP & Canon Cartridges	each	1
All	Std	1900	Bag, Anti-Static 3" X 5" ID	100/pk	1 pk
All	Std	4753	Label Red Adhesive Pull Strip 3/8"X1-1/2", Strip of 10	10/pk	1 pk
PGI-220 Bk	Std	4761	Label Canon PGI-220 Black	each	100
CLI-221 Bk	Std	4764	Label Canon CLI-221 Black	each	100
CLI-221 C	Std	4766	Label Canon CLI-221 Cyan	each	100
CLI-221 M	Std	4767	Label Canon CLI-221 Magenta	each	100
CLI-221 Y	Std	4768	Label Canon CLI-221 Yellow	each	100
CLI-221 Gy	Std	4765	Label Canon CLI-221 Grey	each	100
PGI-225 Bk	Std	4762	Label Canon PGI-225 Black	each	100
CLI-226 Bk	Std	4769	Label Canon CLI-226 Black	each	100
CLI-226 C	Std	4770	Label Canon CLI-226 Cyan	each	100
CLI-226 M	Std	4771	Label Canon CLI-226 Magenta	each	100

CLI-226 Y	Std	4772	Label Canon CLI-226 Yellow	each	100
All	Std	4942	Box White HP50/90/26/29/BC02	each	1
All	Std	LAB104	LAB104, Label "Quality Reman Cart"	each	25
<b>Equipment</b>					
All	Std	1120	EasyFill-M1 *(see video)	each	1
All	Std	6010	Vacuum Generator *(see video)	each	1
All	Std	7075	Scale Digital Compact	each	1

**Instructions:**

\* Go to our website ([www.rjettek.com](http://www.rjettek.com)). Enter the equipment part number in the search field. Click on "Product Description". Scroll down the page to "Videos" and click on the "▶" symbol in the center to start the video.